



# Artificial Urethral Sphincter Dimensions Chart



Model Number	Lumen Diameter	Cuff Width*
AUS-4x11	4mm	11mm
AUS-6x11	6mm	11mm
AUS-8x11	8mm	11mm
AUS-10x11	10mm	11mm
AUS-12x11	12mm	11mm
AUS-14x11	14mm	11mm
AUS-16x11	16mm	11mm
AUS-4x14	4mm	14mm
AUS-6x14	6mm	14mm
AUS-8x14	8mm	14mm
AUS-10x14	10mm	14mm
AUS-12x14	12mm	14mm
AUS-14x14	14mm	14mm
AUS-16x14	16mm	14mm

\*approximate



## The AUS Cuff

To correctly choose the most appropriate size AUS occluder, the diameter of the urethra should be measured by ultrasound. The AUS occluder should fit slightly loose at first; for example if the urethra has an outer diameter of 7mm you would select an occluder with an 8mm lumen diameter.

FELINE PATIENTS- cuff width 11mm  
lumen diameters of 4, 6, 8, 10, 12, 14, & 16mm

CANINE PATIENTS- cuff width 14mm  
lumen diameters of 4, 6, 8, 10, 12, 14, & 16mm

## The Actuating Tube Length

Our stock length is 40cm.

## The Actuating Tube Diameter

- ID .062" x OD .125"