


# Hepatic Shunt Occluders Dimensions Chart



Model Number	Lumen Diameter	Cuff Width*	Cuff Thickness*
HSO-4	4mm	7mm	3mm
HSO-5	5mm	7mm	3mm
HSO-6	6mm	7mm	3mm
HSO-8	8mm	7mm	3mm
HSO-10	10mm	7mm	3mm
HSO-12	12mm	7mm	3mm
HSO-14	14mm	7mm	3mm
HSO-16	16mm	7mm	3mm

*\*approximate*

## The HSO Cuff

The occluder cuff dimensions when deflated are most often the determining factor in the selection of the proper device.

## The Actuating Tube Length

Our stock length is 40cm.

## The Actuating Tube Diameter

- ID .062" x OD .125"

Note: the actuating tube may also be positioned to protrude from the side of the occluder diaphragm (instead of the back) if this alteration is better suited to your application requirements.